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Use

The product is designed for the manufacture of medical and pharmaceutical packaging. The materials comply with the applicable legal and normative requirements.

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Typical Application

Gas-permeable, microbial barrier (sealing component) for medical technology packaging.

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Functional Characteristics
Extremely strong material
Breathable
High tear resistance
Excellent microbial barrier properties

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Storage Condition

Storage in original packaging in a clean, dry environment away from direct exposition to light. Avoid any direct heat. Not to be stored together with solvents or chemicals in the same room.

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Recommended Storage Conditions	
temperature	7°C – 29°C
relative humidity	30%rH-60%rH

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Shelf – life

If stored correctly in accordance with the conditions in this data sheet, the material has a shelf life of 5 years from the date of manufacture.

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If the product is used as part of a sterile barrier system, the user himself is responsible to demonstrate shelf life of the packaging system according to DIN EN ISO 11607-1:2020 for his specific purpose through a validation study.

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Processing Conditions

This material is intended and suitable as a packaging material for medical devices and/or pharmaceutical products. For this reason, the processing conditions for the material are based on the standard environmental conditions for such products (cleanliness, temperature, air humidity) in order to maintain this suitability. In order to determine and specify the effective manufacturing parameters for shaping and/or sealing, the manufacturer is required to carry out a process validation.

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The material is suitable for sealing with films such as PETG, PE, some PP types.

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Biocompatibility

A cytotoxicity test (acc. ISO 10993-5) was carried out for this packaging material and no cytotoxicity was detected.

This sealing component is gas-permeable, but guarantees excellent bacterial and germ resistance (LRV min. 5.2 acc. ASTM F1608), therefore it is suitable for gas sterilization.

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Regulatory Status

MEDIPACK AG does not see itself in the role of a distributor / legal manufacturer within the meaning of the

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Medical Device Regulation EC 2017/745 MDR, but in the role of a manufacturer of packaging components and services. Accordingly, the distributor / legal manufacturer in its role in relation to the MDR is obliged to provide the corresponding tests and evidence for the product and its notified body itself. MEDIPACK AG supports its customers in this respect as far as possible, but the final responsibility remains with the customer - legal manufacturer.

As a Swiss-based manufacturer of medical packaging components and services, we are not directly affected by Regulation 2017/821/EU. Nevertheless, we attach great importance to compliance with this regulation by our raw material suppliers. We require our suppliers to provide the relevant evidence and ensure that their supply chains are free of conflict minerals in order to guarantee the conformity of our entire supply chain.

Conformity

Conformity with DIN EN ISO 11607 - “Packaging for terminally sterilized medical devices”. This product is only one component of the complete sterile barrier system. The conformity of the complete packaging system is the responsibility of the medical device manufacturer.

The product fulfills the basic requirements for material specifications according to DIN EN ISO 11607-1:2020 (Article 5) and is therefore suitable to ensure the sterility and safety of the packaging until the time of use. This applies if the distributor has carried out a process validation in accordance with DIN EN ISO 11607-2:2020.

The product complies with European Directive 94/62/EC (packaging and packaging waste) and Regulation 1935/2004/EC (food contact materials and articles) and its amendments. The materials and raw materials fulfill the requirements of Regulation 10/2011/EU and 2020/1245/EU including all amendments.

The mentioned product contains raw materials that could be produced using substances of animal origin. However, such substances of animal origin are always chemical reaction products in which the animal material is only the original starting material. Animal material as such is not directly present in the product. The relevant regulations of the European Community for the prevention of TSE/BSE contamination have been taken into account and complied with.

With reference to the “Guidance document on minimizing the risk of transmitting animal spongiform encephalopathy agents via human and veterinary medicinal products (EMA/410/01)“ and EC directives, the present material complies with the aforementioned requirements that the raw material undergoes treatments that

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Die Firma MEDIPACK AG sieht sich im Sinne der Medizinprodukte Verordnung EG 2017/745 MDR nicht in der