



STANDARD PEELABLE POUCH FOR HIGH-QUALITY IMPLANTS AND INSTRUMENTS

MEDPOUCH PEEL PA/PE 90

MEDPOUCH PEEL TYVEK® 1073B - PET/PE

validated to ISO 11607


MEDIPACK AG
M E D I Z I N A L V E R P A C K U N G E N

in mm			Art. no. PA/PE	Art. no. Tyvek®- PET/PE
Width (A)	Length (B)	Head (C)		
45	130	15	100.090.01	100.120.01
45	150	15	100.090.02	100.120.02
45	240	15	100.090.03	100.120.03
45	340	15	100.090.04	100.120.04
45	490	15	100.090.05	100.120.05
45	600	15	100.090.06	100.120.06
55	150	15	100.091.01	100.121.01
65	130	25	100.092.01	100.122.01
65	150	25	100.092.02	100.122.02
65	194	25	100.092.07	100.122.07
65	240	25	100.092.03	100.122.03
65	340	25	100.092.04	100.122.04
65	490	25	100.092.05	100.122.05
65	570	25	100.092.06	100.122.06
65	600	25	100.092.09	100.122.09
65	650	25	100.092.08	100.122.08
85	150	30	100.093.01	100.123.01
85	180	30	100.093.02	100.123.02
85	220	30	100.093.03	100.123.03
85	240	30	100.093.04	100.123.04
85	340	30	100.093.05	100.123.05
85	490	30	100.093.06	100.123.06
85	650	30	100.093.07	100.123.07
100	200	35	100.110.02	100.130.01
100	280	35	100.110.03	100.130.02
100	650	35	100.110.01	100.130.03
100	720	35	100.110.04	100.130.04
110	205	35	100.094.01	100.124.01
110	250	35	100.094.04	100.124.04
110	340	35	100.094.02	100.124.02
110	380	35	100.094.03	100.124.03
120	220	35	100.095.01	100.125.01
125	650	35	100.111.01	100.131.01
130	220	35	100.096.01	100.126.01
130	340	35	100.096.02	100.126.02
150	220	40	100.097.01	100.127.01
150	265	40	100.097.02	100.127.02
150	340	40	100.097.03	100.127.03
150	430	40	100.097.04	100.127.04
150	450	40	100.097.05	100.127.05
150	490	40	100.097.06	100.127.06
175	240	40	100.098.01	100.128.01
175	285	40	100.098.02	100.128.02
175	340	40	100.098.03	100.128.03
175	380	40	100.098.04	100.128.04
200	380	40	100.099.01	100.129.01
200	490	40	100.099.02	100.129.02

STANDARD PEELABLE POUCH

Standard MEDPOUCH peeling pouches are suitable for the packaging of high-quality implants and instruments. They can be sterilised by gamma irradiation or ethylene oxide gas (Tyvek® pouches). Pouch production takes place under clean room conditions in accordance with ISO Class 7.


Quality MEDPOUCH PEEL PA/PE 90

Frontside: PA/PE 30/70 with PA/PE 30/60 peel

Quality MEDPOUCH TYVEK® 1073B – PET/PE

Frontside: PET/PE-film 12/50my

Backside: Tyvek® 1073B 63my uncoated



More peelable pouches are available by request.

MEDIPACK AG

Mühlentalstrasse 184-188
8200 Schaffhausen
Schweiz
Telefon: +41 52 630 36 36

info@medipack.ch

MEDIPACK AG

Zweigniederlassung Deutschland
c/o Dachser SE
Josef-Schüttler-Str. 5
78224 Singen
Deutschland

info@medipack.ch



Certified to ISO 13485:2016
Technical subject to change

www.medipack.ch